

OLYMPUS[®]

RHINO-LARYNGO FIBERSCOPE
OLYMPUS ENF TYPE XP




The Optimal Choice for Observation of Narrow Upper Airway Passages




Ø1.8mm

Ultra-Slim 2.2 mm Insertion Tube—Perfectly Sized for Narrow Nasal Cavities.

Featuring an ultra-slim insertion tube diameter of 2.2 mm as well as high resolution optics and brightness, the new Olympus ENF TYPE XP is perfectly suited for observing the nasal cavity, larynx and pharynx. Optimally sized for narrow nasal cavities, it also makes observation of adult paranasal sinuses easier than ever, while its angulation range of 130° both up and down enables trachea insertion maneuverability. For minimally invasive insertion and clear, reliable observation, depend on the Olympus ENF TYPE XP Rhino-Laryngofiberscope.



Ultra-Slim 2.2 mm Insertion Tube Diameter
Optimally sized for smaller nasal cavities. This endoscope also makes observation of adult paranasal sinuses easier than ever.




Wide Angulation Range—130° Up and Down
The endoscope wide bending range of 130° up/130° down allows trouble-free insertion, even through the trachea, allowing you to focus more on observation.

Specifications

Field of view	75°
Depth of field	2.5~50 mm
Range of tip bending	Up 130°/Down 130°
Distal end outer diameter	ø1.8 mm
Insertion tube diameter	ø2.2 mm
Working length	300 mm
Total length	530 mm

Olympus Rhino-Laryngofiberscope Line-Up

Versatile, durable and reliable Olympus rhino-laryngofiberscopes ensure smooth operation. A full range of rhino-laryngo fiberscopes for your procedural needs. Olympus rhino-laryngofiberscope—the right choice for superior performance.



ENF-P4



ENF-T3



ENF-L3

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS OPTICAL CO., LTD.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
OLYMPUS WINTER & IBE GMBH
Kuehnstrasse 61, 22045 Hamburg, Germany
OLYMPUS OPTICAL CO. (EUROPA) GMBH.
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
OLYMPUS LATIN AMERICA, INC.
6100 Blue Lagoon Drive, Suite 390, Miami, Florida 33126-2087, U.S.A.
OLYMPUS KEYMED
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
491B, River Valley Road, #12-01/04, Valley Point Office Tower, Singapore 248373
OLYMPUS HONG KONG AND CHINA LIMITED.
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS BEIJING REPRESENTATIVE OFFICE
Room 1401, East Tower, Gong Yuan No.6 Royal Palace, No.6 Gong Yuan Xijie, JianGuaMenNei, DongCheng District, Beijing, 100005, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl. 2, Russia
OLYMPUS AUSTRALIA PTY. LTD.
31 Gilby Road, Mount Waverley, VIC., 3149, Australia