

OLYMPUS

Your Vision, Our Future


RIGID LARYNGOSCOPY

HD Laryngoscope


ENT

HD²

High Definition · High Durability



HD² LARYNX



HIGH DEFINITION




High image quality – great convenience

Be on the safe side:
The HD laryngoscope offers a unique HD optical system. Contrast and colour reproduction are at the highest level of technology. Finest structures visible!


Enjoy the great convenience:
Thanks to the high depth of field, no refocusing is required during procedures.

HDTV¹⁰⁸⁰ⁱ
EVIS
EXERA II

HYPERGLIDE LARYNGOSCOPY



HIGH DURABILITY



Innovative design – great durability

Reduce costs over lifetime:
The new patented lens system ensures serious durability.

Shocks are absorbed by the innovative lens protection system.

Lens breakage is reduced during daily work and reprocessing.



Designed for patients' convenience

Less pharyngeal reflex:
The distance between the optics and the distal end is reduced to a minimum.



Ready for use

Whatever is required:
The HD laryngoscope is compatible with relevant reprocessing techniques – autoclaving, automated cleaning and disinfection.


Detach the handle for easy and complete cleaning.

HD² LARYNGOSCOPY


High Definition Laryngoscope

HD Laryngoscopes are delivered with detachable handle.

WA96100A Telescope, HD,
10 mm,
167 mm working length,
70° direction of view,
autoclavable




WA96105A Telescope, HD,
10 mm,
167 mm working length,
90° direction of view,
autoclavable



Accessories


WA03210A Light-guide cable,
size M,
plug type,
3 m,
CF type



A3201 Adapter,
Olympus light-guide cable to light sources
Storz,
Aesculap,
Linvatec




WA05970A Instrument tray with lid,
with silicone mat WA05980A,
outer dimensions: 537 x 139 x 268 mm



Video System


Video system center "CV-180"



Camera head,
HDTV 1080i, angled,
with three remote
control switches,
autoclavable,
with integrated video adapter
(clamping connector),
NBI compatible



N3046640 "OTV-S7ProH-HD-10E", 1.0 x
N2638640 "OTV-S7ProH-HD-12E", 1.2 x



OLYMPUS

Manufactured by
OLYMPUS WINTER & IBE GMBH

Kuehnstraße 61, 22045 Hamburg, Germany

Distributed by
OLYMPUS MEDICAL SYSTEMS CORPORATION
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
OLYMPUS SURGICAL & INDUSTRIAL AMERICA INC.
KEYMED (MEDICAL & INDUSTRIAL EQUIPMENT) LTD.
OLYMPUS SINGAPORE PTE LTD.
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
OLYMPUS AUSTRALIA PTY. LTD.
OLYMPUS LATIN AMERICA INC.
OLYMPUS KOREA CO., LTD.
OLYMPUS (BEIJING) SALES & SERVICE, CO., LTD.

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
Wendenstraße 14-18, 20097 Hamburg, Germany
One Corporate Drive, Orangeburg, New York 10962, U.S.A.
KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
491B, River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
31 Gilby Road, Mount Waverley, Victoria 3149, Australia
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126-2097, U.S.A.
4F, Gyeongam Bldg., 157-27, Samseong-Dong, Kangnam-Gu, Seoul 135-090, Korea
R1202, NCI Tower, A12 Jianguomenwai Dajie, Chaoyangqu Beijing 100022, China