

PRECISION TOOLS FOR LASER TREATMENT OF BPH



SIDE-FIRE LASER




The 22.5 Fr. OES Pro continuous flow laser-cystoscope was specifically designed for side-fire laser treatment of the prostate – ideal for both Laserscope (GreenLight PV System)* or Lumenis/Boston Scientific (VersaPulse DuoTome)**.



Interstitial laser treatment

- Instrument guide right to the tip
- Controlled flow characteristics
- Alternative use: collagen injection



OES 4000 Continuous Flow Cystoscope

ER TREATMENT



- Specially designed working insert WA22820A for GreenLight Photo-selective Vaporization of the Prostate (PVP) maximises control of the laser fiber due to the separate guidance channel
- High-flow working insert WA22850A for both, PVP procedures or holmium laser ablation of the prostate

Unique laser control adapter for GreenLight fiber ensures best longitudinal control and minimises sheath or scope tip damage.

* Laserscope and GreenLight PV are registered trademarks of Laserscope
** Lumenis, VersaPulse and DuoTome are registered trademarks of the Lumenis group of companies


Alternative Option: Side-Fire treatment with resection sheath

- 26 Fr. continuous flow sheath
- 8 Fr. channel


OES Pro System




Side-Fire Laser Cystoscope




Side-Fire Laser-Treatment



Interstitial Laser Cystoscope



HoLEP Laser Resectoscope



OLYMPUS

Manufactured by
OLYMPUS WINTER & IBE GMBH

Kuehnstraße 61, 22045 Hamburg, Germany

Distributed by
OLYMPUS MEDICAL SYSTEMS CORPORATION
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
OLYMPUS SURGICAL & INDUSTRIAL AMERICA INC.
OLYMPUS AMERICA INC.
KEYMED (MEDICAL & INDUSTRIAL EQUIPMENT) LTD.
OLYMPUS SINGAPORE PTE LTD.
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
OLYMPUS AUSTRALIA PTY. LTD.
OLYMPUS LATIN AMERICA INC.
OLYMPUS KOREA CO., LTD.
OLYMPUS (BEIJING) SALES & SERVICE, CO., LTD.

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
Wendenstraße 14-18, 20097 Hamburg, Germany
One Corporate Drive, Orangeburg, New York 10962, U.S.A.
Two Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
491B, River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
31 Gilby Road, Mount Waverley, Victoria 3149, Australia
6100 Blue Lagoon Drive, Suite 390 Miami, FL 33126-2087, U.S.A.
8F, Hyundai Marines Bldg., 646-1, Yeoksam-Dong, Kangnam-Gu, Seoul 135-080, Korea
R1202, NCI Tower, A12 Jianguomenwai Dajie, Chaoyangqu Beijing 100022, China