

**OLYMPUS**


Your Vision, Our Future

UROLOGY

**Long Resectoscope**

OES Pro


# THE EXTRA LENGTH TO MEET ALL YOUR NEEDS



**7 cm more**  
**2.7 inches more**

# OLYMPUS LONG OES PRO RESECTOSCOPE – THE SHORT SOLUTION TO A LONG PROBLEM

The long resectoscope was developed for the resection of bladder tumors as well as for interventions in patients with long urethras and semi-flexible penile implants.




- ✓ The modular design of the long resectoscope allows for 24 Fr. single flow, or 27 Fr. continuous flow use.
- ✓ Two different, specially adapted outer sheaths allow for tumor resection in the bladder and for resection in the prostate.
- ✓ The long resectoscope with its continuous flow system provides improved visualization under low pressure conditions.

THE LONG RESECTOSCOPE – THE RIGHT FIT

# OES PRO LONG RESECTOSCOPE

## Long Telescopes

WA20021A Telescope, 4 mm, long, 12° direction of view, autoclavable




WA20023A Telescope, 4 mm, long, 70° direction of view, autoclavable




To be used in connection with bridge A22093A

## Working Elements

WA22066A Working element, active




WA22067A Working element, passive




A0393 HF cable, monopolar, 3.5 m, for UES-30/-40, Erbe ICC Int., Valleylab (new)  
 A0391 for UES-10/-20 and Valleylab (old)  
 A0392 for Martin and Berchtold

## HF-Resection Electrodes


WA22037A HF-resection electrode, loop, 12°, reusable, 1 pc.  
 WA22037C sterile, single use, 12 pcs.



WA22038A HF-resection electrode, roller, 12°, reusable, 1 pc.  
 WA22038C sterile, single use, 12 pcs.




WA22039A HF-resection electrode, angled loop, 12°, reusable, 1 pc.  
 WA22039C sterile, single use, 12 pcs.




## 24 Fr. Standard/Inner Sheaths

WA22017A Inner sheath, 24 Fr., incl. standard obturator



WA22017T Inner sheath, 24 Fr., incl. obturator with deflecting tip



Irrigation port, for resection sheaths,  
 A22051A 1 stopcock, rotatable  
 A22052A 1 Luer-lock connector, rotatable  
 A22053A 2 stopcocks, rotatable  
 A22054A 1 vertical stopcock, fixed  
 A22055A 1 vertical Luer-lock connector, fixed




WA22016A Obturator, optical, for telescope WA20021A




## 27 Fr. Outer Sheaths

WA22018A Outer sheath, 27 Fr., for the resection of bladder tumors, working length: 262 mm




WA22019A Outer sheath, 27 Fr., for the resection of prostatic tissue, working length: 262 mm



## Bridge

A22093A Bridge



For inspection of the bladder with telescope WA20023A

## OES Pro Light-Guide Cable


WA03200A Light-guide cable, size S, plug type, 3 m



A0460 Light-guide adapter, for OES 4000 light-guide cable  
 O0332 for Wolf light-guide cable



## System Chart



**OLYMPUS**

Manufactured by  
**OLYMPUS WINTER & IBE GMBH**

Distributed by  
**OLYMPUS MEDICAL SYSTEMS CORPORATION**  
**OLYMPUS MEDICAL SYSTEMS EUROPA GMBH**  
**OLYMPUS SURGICAL & INDUSTRIAL AMERICA INC.**  
**OLYMPUS AMERICA INC.**  
**KEYMED (MEDICAL & INDUSTRIAL EQUIPMENT) LTD.**  
**OLYMPUS SINGAPORE PTE LTD.**  
**OLYMPUS MOSCOW LIMITED LIABILITY COMPANY**  
**OLYMPUS AUSTRALIA PTY. LTD.**  
**OLYMPUS LATIN AMERICA INC.**  
**OLYMPUS KOREA CO., LTD.**  
**OLYMPUS (BEIJING) SALES & SERVICE, CO., LTD.**

Kuehnstraße 61, 22045 Hamburg, Germany

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan  
 Wendenstraße 14-18, 20097 Hamburg, Germany  
 One Corporate Drive, Orangeburg, New York 10962, U.S.A.  
 Two Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.  
 KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom  
 491B, River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore  
 117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia  
 31 Gilby Road, Mount Waverley, Victoria 3149, Australia  
 6100 Blue Lagoon Drive, Suite 390 Miami, FL 33126-2087, U.S.A.  
 8F, Hyundai Marines Bldg., 646-1, Yeoksam-Dong, Kangnam-Gu, Seoul 135-080, Korea  
 13F - 14F, East Tower, Gongyuan No. 6 Royal Palace, No. 6 Gongyuan Xijie, Jianguomennei, Dongcheng District, Beijing 100005, China