


Variety in Mini-Hysteroscopes



1.9 MM MINI SYSTEM OES 4000 HYSTEROSCOPES



Focussing on minimal invasiveness only flexible hysteroscopes could meet the expectations for maximal patient friendliness. The new Olympus 1.9 mm system combines the benefits

of flexible as well as rigid hysteroscopes: patient comfort (slim diameter), easy reprocessing (autoclavability) and efficient therapy (5 or 7 Fr. working channel).



Single flow:
3 mm sheath for best patient comfort



Durability by semiflexible optical fibers



Continuous flow:
4.5 mm sheath including 5 Fr. working channel



All-in-one system for diagnosis and therapy








Continuous flow:
4.5 mm and 7 Fr. working channel



No anesthesia required

Technical Data

A4677A A4676A	Telescope, 1.9 mm, autoclavable, 30° direction of view 0° direction of view	
A4775	Sheath, single flow, 3 mm	
A4776	Sheath, continuous flow, 4.5 mm, 5 Fr. channel	
A46700A	Compact hysteroscope, 4.5 mm, 12°, continuous flow, 7 Fr. channel	
WA03200A	Light-guide cable, 3.0 mm, plug type, CF type	

OLYMPUS

Manufactured by
OLYMPUS WINTER & IBE GMBH

Distributed by
OLYMPUS MEDICAL SYSTEMS CORPORATION
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
OLYMPUS SURGICAL & INDUSTRIAL AMERICA INC.
OLYMPUS AMERICA INC.
OLYMPUS KEYMED GROUP LTD.
OLYMPUS SINGAPORE PTE LTD.
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
OLYMPUS AUSTRALIA PTY. LTD.
OLYMPUS LATIN AMERICA INC.
OLYMPUS KOREA CO., LTD.
OLYMPUS (BEIJING) SALES & SERVICE, CO., LTD.

Kuehnstraße 61, 22045 Hamburg, Germany

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
Wendenstraße 14-18, 20097 Hamburg, Germany
One Corporate Drive, Orangeburg, New York 10962, U.S.A.
Two Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
491B, River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
31 Gilby Road, Mount Waverley, Victoria 3149, Australia
6100 Blue Lagoon Drive, Suite 390 Miami, FL 33126-2087, U.S.A.
8F, Hyundai Marines Bldg., 646-1, Yeoksam-Dong, Kangnam-Gu, Seoul 135-080, Korea
13F - 14F, East Tower, Gongyuan No. 6 Royal Palace, No. 6 Gongyuan Xijie, Jianguomennei, Dongcheng District, Beijing 100005, China